


720 Second Ave, Eau Claire, WI 54703 www.echealthdepartment.org P: 715-839-4718 F: 715-839-1674

## **Thermometer Calibration Log**

Calibrate all thermometers at least once a week and when dropped



## **How to Calibrate a Thermometer**

- 1. Fill a glass with ice (crushed is best) and add a little water.
- 2. Submerge thermometer stem.
- 3. Wait a minute for temperature to adjust.
- 4. If thermometer does not read 32°F when in the ice water, calibrate:
  - a. If dial-use a wrench to turn the nut under the dial face until thermometer reads 32°F.
  - b. If digital-see manufacturer instructions or send it to the manufacturer to calibrate.

Note: Keep the thermometer in the water when calibrating.

	Enter initials when thermometer is calibrated correctly				
Date	Thermometer 1	Thermometer 2	Thermometer 3	Thermometer 4	Thermometer 5

Replace thermometers if they cannot be accurately calibrated

Reviewed by	
-------------	--